

L8 ANSWER 4 OF 9 CAPLUS COPYRIGHT 2000 ACS
TI Ionization and tautomerism of **fluorescein**, **rhodamine**
B, N,N-diethylrhodol and related dyes in mixed and nonaqueous solvents
AU Mchedlov-Petrussyan, Nikolay O.; Kulhtik, Valentina I.; Alekseeva, Vera I.
SO Dyes Pigm. (1994), 24(1), 11-15
CODEN: DYPYDH; ISSN: 0143-7208
AB The prot. equiv. of **fluorescein**, **rhodamine** B and
of the aryl amino-oxyanthene dye, N,N-diethylrhodol (a 'hybrid' of
rhodamine B and **fluorescein**) were studied in aq. DMSO
and EtOH (9:1 org. cosolvent). The pKa values of these dyes, as
well as of related substances, were detd. On the basis of the visible
absorption spectra in various solvents conclusions were made about
tautomerism in the dye mols. Values of the tautomeric equil. consts. and
of the microscopic ionization consts. were obtained. Some new data on the
tautomerism of oxyanthene monoanions in MeOH were presented.

=> d t i a s c a i s s

L8 ANSWER 3 OF 9 MEDLINE
TI Phagosomal pH determination by dual fluorescence flow cytometry.
AU Vergne I; Constant D; Lameille S
SO ANALYTICAL BIOCHEMISTRY, (1998 Jan 1) 255 (1) 127-32.
Journal Code: 4NK. ISSN: 0003-2697.
AB Several methods have been developed to measure the pH of phagosomes, using
fluorescein derivatives as reporter of pH, and spectrofluorimetry,
fluorescence microscopy, or flow cytometry as quantification technique.
All have major disadvantages, including either a slow or inaccurate
response. In the present study, pH determination was achieved on J774-cell
phagosomes containing dual-labeled zymozan particles using dual
fluorescence flow cytometry with an argon-ion laser excitation wavelength
at 488 nm. This allowed zymozan-containing macrophages to be distinguished
from other cells and their fluorescence to be measured rapidly. The use of
a new probe, namely Oregon Green 488 which has a pKa lower than
carboxyfluorescein with the same maximum excitation and emission
wavelengths, allowed investigation of pH value below 5. The dual labeling
with Oregon Green 488 and carboxytetramethylrhodamine as pH-sensitive and
pH-insensitive probes, respectively, overcame the absence of an isobestic
point in the Oregon Green 488 spectrum. The phagosomal pH was determined
using a calibration curve of phagosomal pH established by adding
ionophores to phagocyte suspension and measuring the fluorescence
intensity ratio (535 nm/585 nm) for different pHs. A phagosomal pH of 4.5
± 0.1 can be accurately determined. This method permits pH measurements
down to 4, even in the presence of nonengulfed reporter particles.

JP 10.1.98